

Forano, Front Part, DIN, Ar, Colour-Neutral



Manufacturer: Greggersen

Product 905117

Number:

Description

Terminal Unit FORANO - Front Part DIN

Colour-Neutral - Version III from 08/2019

Gas type: Argon / Ar

Labelling: argon/ Ar

Gas-specific Socket Assembly for installation
on Base Part FORANO.

Socket Assembly according to DIN 13260 Part 2
with parking and operating position.

Release via Release Ring.

Colour coding:

Colour-Neutral (black)

Release Ring labelling: printed

Language: English

Cartridge incl. sealing elements

Technical data:

Gas type: Argon/ Ar

Labelling: argon/ Ar

Colour coding: Colour-Neutral (black)

Front diameter: Ø 42mm

Scope of delivery:

- gas-specific Socket Assembly
- Release Ring incl. print labelling
- Cartridge with sealing elements

Properties

Gas Type

Ar-Argon
