

Falke without Flow easy O2 NF 3/8"



Manufacturer: Greggersen

Product 910389

Number:

Description

Falke "easy", 450 kPa, G3/8",
Compact Pressure Regulator with Cylinder Connection
according to DIN 477, contains gauge, free flow rate at a
pressure of 450 kPa.

Gas type: Oxygen / O2

Connection: Cylinder Connection according to NF N29-650, F-627

Outlet: G3/8" with Non-Return Valve, bottom

Output: 450 kPa

Properties

Gas Type

O2-Oxygen