

Hose, VAC, 5m, DIN - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 900473

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe,

Hose and Device Connection NIST

Specification:

Gas type: Vacuum / VAC

Inlet: Probe according to DIN 13260 Part 2

Hose length: 5.0 m

Colour coding: Colour-Neutral

Outlet: gas-specific screw connection NIST

Properties

Gas Type

VAC-Vacuum