

Varioport Pirol 2.0 -90, VAC, NIST



Manufacturer:	Greggersen
----------------------	------------

Product	904609
----------------	--------

Number:	
----------------	--

Description

Portable suction unit Varioport with Vacuum Regulator Pirol -90

Medical, portable suction unit for vacuum
to accommodate 2 secretion jars.

The holder for the jars is universal. Standard jars in
sizes 0.5 to 3.0 L from various manufacturers can be used
(e.g. Serres, Abbott, Medela, Ardo, Cardinal, and many others).

The carrying frame features a sturdy handle and a mount for
a 25x10 mm equipment rail. The vacuum control unit is located
in the front area, allowing the operator quick and direct access.

The gas inlet as a NIST connection is located at the rear of
the unit, so that the supply hose remains outside the working area.

Scope of Delivery:

- Vacuum Regulator Pirol-90 with fine dosing valve, rotatable Vacuum Gauge
- Carrying frame, coated
- 2 pcs highly flexible retaining collars

Technical Data:

Gas type: Vacuum - VAC

Inlet: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF with vacuum tubing nipple

Control range: 0 to -90 kPa $\pm$ 10%

Flow Rate: > 20 l Freeflow/min

Weight: approx. 3 kg (without jars)

Properties

Gas Type	VAC-Vacuum
-----------------	------------
