


Forano, Front Part, DIN,AIR-800, Colour-Neutral



Manufacturer: Greggersen

Product 905119

Number:

Description

Terminal Unit FORANO - Front Part DIN

Colour-Neutral

Gas type: Compressed Air-800 / AIR-800

Labelling: medical air 800 / AIR-800

Gas-specific Socket Assembly for mounting
on Base Part FORANO.

Socket Assembly according to DIN 13260 Part 2
with parking and operating position.

Release via Release Ring.

Colour coding:

Colour-Neutral (black)

Release Ring labelling: printed

Language: English

Cartridge incl. sealing elements

Technical data:

Gas type: Compressed Air-800 / AIR-800

Labelling: medical air 800 / AIR-800

Colour coding: Colour-Neutral (black)

Front diameter: Ø 42mm

Scope of delivery:

- gas-specific Socket Assembly
- Release Ring incl. printed labelling
- Cartridge with sealing elements

Caution – Compatibility Notice

Before using the AIR-800 Terminal Unit, it must be ensured that the Probes used comply with the requirements of DIN 13260-2:2025.

In the past, similarly designed Probes were distributed that are not compatible with the AIR-800 Terminal Unit. The compatibility of the Probes used must therefore be checked prior to use.

Properties

Gas Type

AIR-800-Compressed Air 8bar
