

Vacuum Regulator Spatz-90 Probe - BS



Manufacturer: Greggersen

Product 901509

Number:

Description

Vacuum Regulator Spatz-90, Plug-in Unit BS

Medical vacuum regulator based on the bypass principle,
with vacuum gauge, fine adjustment and quick-action shut-off valve.

Technical Data:

Gas type: Vacuum / VAC

Inlet: Probe to BS 5682:1998

Control range: -20 to -85 kPa (-0.2 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple