


Falke without Flow N2O PIN +1DIN



Manufacturer:	Greggersen
Product	500030
Number:	

Description

Falke, without flow with 1 quick-action coupling,
 Compact pressure regulator with PIN INDEX high-pressure connection,
 Contents gauge, with one quick-action coupling at a
 pressure of 450 kPa.

Gas type: Nitrous oxide / N2O

Connection: PIN INDEX high-pressure connection according to CGA V-1; No:910

Outlet: Quick-action coupling according to DIN 13260 Part 2, top

Output: 450 kPa

Special version

PIN INDEX high-pressure connection

Properties

Gas Type	N2O-Nitrous Oxide
-----------------	-------------------