

Kolibri 2 DIN Probe O₂ Paediatric



Manufacturer: Greggersen

Product 906402

Number:

Description

KOLIBRI 2, O₂, 0...5 l/min

- Paediatric DIN PLUG-IN UNIT -

Flow Settings:

0 / 0.1 / 0.3 / 0.5 / 0.6 / 0.8 / 1 / 2 / 3 / 4 / 5 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selector device

Medical compact flowmeter for connection

to a terminal unit, with handwheel – incrementally

adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and readability

of the selected flow setting from various viewing directions

and angles.

Frontal operating concept -> increased safety when
setting the correct flow rate

Gas Type: Oxygen / O₂

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow Rate: 0...5 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:

900619 - Hose Connection 9/16-18 UNF

Properties

Gas Type	O ₂ -Oxygen
Inlet Standard	DIN
Variant	Paediatric