

Hose, AIR, 1.5m, NF - NIST, ISO



Manufacturer: Greggersen

Product 902514

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow limiter,
hose and NIST equipment connection

Specification:

Gas type: Compressed Air / AIR

Inlet: Probe according to NF S 90-116

Hose length: 1.5 m

Colour coding: ISO 32 (black/white)

Outlet: gas-specific screw connection NIST

Properties

Gas Type AIR Compressed Air