

Vacuum Reg. Skua 2.0 -30,AIR,PROBE,DIN



Manufacturer: Greggersen

Product 904602

Number:

Description

Compressed-Air Vacuum Regulator Skua -30, Plug-in Unit DIN

Medical compressed-air vacuum regulator for dosing vacuum.

The vacuum is generated by means of compressed air according to the Venturi principle.

With 360° rotatable vacuum gauge; large, ergonomic adjustment handwheel; rotatable exhaust sleeve.

The adjustment mechanism enables stepless, precise and intuitive regulation of the vacuum. The rotatable vacuum gauge allows the user to read the set value from an optimal viewing angle. The exhaust sleeve also allows the air flow to be individually directed.

The regulator unit complies with the requirements of DIN EN ISO 10079-3 (high vacuum /high flow).

Material:

Housing: aluminium, anodised

Handwheel: plastic

Probe front part: stainless steel

Technical data:

Gas type: med. compressed air / AIR

Inlet: probe in accordance with DIN 13260 Part 2

Control range: 0 to -30 kPa $\pm$ 10%

Flow rate: >20 l freeflow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------
