

Portable Emergency Supply Habicht Tube



Manufacturer: Greggersen

Product 902696

Number:

Description

Emergency supply unit,
portable, without oxygen cylinder
- with tube flowmeter -

Emergency supply unit, oxygen flowmeter
and ejector for vacuum generation
consisting of:

Stainless steel carrying frame with stand, to accommodate
a 3 l O₂ cylinder, suitable
for mounting on equipment rail or bed,
central supply block with NIST inlet,
consisting of:

Flowmeter with fine dosing valve for regulating the flow
and ejector with vacuum gauge and fine dosing valve;

Gas type: Oxygen / O₂

Connection: NIST in accordance with DIN EN ISO 18082

Outlet: Tubing nipple / vacuum tubing nipple

Performance: FLOW : 0...15 l/min

VAC: adjustable between
0 and -85 kPa (0 to -0.85 bar)