

Kolibri, DIN Probe, O2 10 Flows



Manufacturer: Greggersen

Product 901406

Number:

Description

Kolibri, O2, 10 Flows

-Plug-in Unit-

Medical Compact Flowmeter according to

DIN EN ISO 15002, for coupling

into a terminal unit, with handwheel

stepwise adjustable flow rate,

outlet 9/16-18 UNF and tubing nipple

Gas type: Oxygen

Connection: Probe according to DIN 13260 Part 2

Total length: 88 mm (probe front part to centre of tubing nipple outlet)

Outlet: 9/16-18 UNF / Tubing Nipple

Flow rate: 0... l/min

Flow settings: +++ Please specify in l/min +++

Properties

Gas Type	O2-Oxygen
----------	-----------
