

Vacuum Regulator Spatz -16 DIN Probe long



Manufacturer: Greggersen

Product 900898

Number:

Description

Vacuum Regulator Spatz-16, Plug-in Unit DIN
extended version (+57 mm)

Medical vacuum regulator with vacuum gauge,
fine dosing and quick-shut-off valve.

Technical data:

Gas type: Vacuum / VAC

Inlet: Probe according to DIN 13260 Part 2

Control range: 0 to -16 kPa (0 to -160 mbar)

Outlet: 9/16-18 UNF with vacuum tubing nipple