

Vacuum Regulator Specht -16 AIR DIN Probe



Manufacturer: Greggersen

Product 900958

Number:

Description

Compressed-Air Vacuum Regulator Specht-16, Plug-in Unit DIN

Medical compressed-air vacuum regulator for vacuum generation, with vacuum gauge, fine dosing and quick-action shut-off valve, silencer and filter mount.

Technical data:

Gas type: Compressed air / AIR

Inlet: Probe according to DIN 13260 Part 2

Control range: 0 to -16 kPa (0 to -160 mbar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type AIR Compressed Air