


Kolibri, Rail CO2 Standard - NIST



Manufacturer: Greggersen

Product 902693

Number:

Description

Kolibri, CO₂, 0...15 l/min

RAIL-MOUNTED UNIT NIST

Medical compact flowmeter according to DIN EN ISO 15002;

Flow rate is set via handwheel in 10 flow settings;

Gas inlet NIST integrated into rail clamp;

Outlet 9/16-18 UNF and tubing nipple

Gas type: Carbon dioxide / CO₂

Connection: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF / Tubing Nipple

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

