


Forano, Front Part, DIN, AIR ,ISO



Manufacturer: Greggersen

Product 905101

Number:

Description

Terminal Unit FORANO - Front Part DIN

ISO 32 - Version III from 08/2019

Gas type: Compressed Air / AIR

Labelling: medical air / AIR

Gas-specific Socket Assembly for mounting
on FORANO Base Part.

Socket Assembly according to DIN 13260 Part 2
with parking and operating position.

Release via Release Ring.

Colour coding:

according to DIN EN ISO 9170-1 (ISO32 - black/white)

Release Ring labelling: printed

Language: English

Cartridge incl. sealing elements

Technical Data:

Gas type: Compressed Air / AIR

Labelling: medical air / AIR

Colour coding: according to ISO 9170-1 (black/white)

Front diameter: Ø 42mm

Scope of Delivery:

- Gas-specific Socket Assembly
- Release Ring incl. printed labelling
- Cartridge with sealing elements

Properties

Gas Type	AIR Compressed Air
----------	--------------------
