

Forano, Base Part, Straight Pipe, DIN, VAC



Manufacturer: Greggersen

Product 903202

Number:

Description

Terminal Unit Forano Base Part DIN/AGA/NF/BS

Gas-specific base part for installation in
supply rails or ducts.

Base part with spacer washer and mounting strip;
brazed copper pipe, integrated non-return and maintenance valve.

Scope of delivery:

- gas-specific base part with Cu pipe
- spacer washer
- mounting strip

Technical data:

Gas type: Vacuum / VAC

Standard AGA: UTL (AGSS)

Installation type: Media supply rail

Inlet: Cu pipe 8x1 mm straight, top

Properties

Gas Type VAC-Vacuum