

Kolibri 2 Rail CH NIST AIR Standard



Manufacturer: Greggersen

Product 906446

Number:

Description

KOLIBRI 2, AIR, 0...15 l/min

- Standard CH RAIL-MOUNTED UNIT -

Flow settings:

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selector device

Medical compact flowmeter for mounting on an

equipment rail; gas inlet NIST integrated in the rail clamp;

with handwheel - incrementally adjustable flow rate, outlet 9/16-18 UNF thread

-> With special, spring-loaded

rail clamp for rails up to 50mm <-

Ergonomically optimised adjustment device and readability of the set flow setting from various viewing directions and angles possible.

Frontal operating concept -> increased safety when setting the correct flow rate

- Gas type: Med. compressed air / AIR
- Connection: NIST according to DIN EN ISO 18082
- Outlet: 9/16-18 UNF
- Flow rate: 0...15 l/min at 450 kPa (4.5 bar) inlet pressure

Accessories required for operation:
900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	AIR Compressed Air
Inlet Standard	CH NIST
Variant	Standard