


Falke with Flow Standard O2 DIN 9/16 AL long



Manufacturer: Greggersen

Product 505100

Number:

Description

Falke with Flow, Standard,
Compact pressure regulator with contents gauge;
step-adjustable flow rate

Gas type: Oxygen / O2

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 9

Outlet: 9/16-18 UNF (tubing nipple), bottom

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

STEM LENGTH 56 mm

Properties

Gas Type O2-Oxygen