

Hose, AIR, 1.5m, MEDAP - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 500154

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe,

hose and equipment connection NIST

Specification:

Gas type: Compressed Air / AIR

Inlet: MEDAP

Hose length: 1.5 m

Colour coding: Colour-Neutral

Outlet: gas-specific screw connection NIST

Properties

Gas Type

AIR Compressed Air