

## Flowmeter, DIN, O<sub>2</sub>, 0-15l/min antimagnetic



**Manufacturer:** Greggersen

**Product** 900772

**Number:**

### Description

Flowmeter, O<sub>2</sub>, 0...15 l/min, non-magnetic

-PLUG-IN UNIT-

Medical flowmeter, pressure-compensated according to  
DIN EN ISO 15002, with fine adjustment valve for  
dosing the flow, for direct connection to  
a terminal unit

Version: non-magnetic

\*\*\*\*\*

Gas type: Oxygen - O<sub>2</sub>

Connection: Probe according to DIN 13260 Part 2

Overall length: 86.5 mm (probe front part to centre of flow outlet)

Outlet: 9/16 - 18 UNF

Flow rate: 0...15 l/min continuously adjustable

Inlet pressure: 450 kPa +/- 50 kPa

Properties

|          |           |
|----------|-----------|
| Gas Type | O2-Oxygen |
|----------|-----------|

---