

Automatic Manifold MediControl MC2050R



Manufacturer: Greggersen

Product 326052

Number:

Description

Automatic Manifold MediControl 2050R

Flow rate: 50 m³/h

Inlet pressure: max. 200 bar

Outlet pressure: 5 bar +/- 0.5 bar

Electr. connection: 230 V AC

System language: German

Mode: 2 cylinder banks and reserve

Gas type: ++++++

Consisting of:

1 x MC2050

1 x Pipe Set MC2050 + Reserve Inlet Pipe

Microprocessor-controlled gas supply system for 2-sided cylinder supply, design according to DIN EN ISO 7396-1, LCD display with indication of operating status and service messages in plain text, forwarding of

messages via potential-free contacts, sensor-based pressure monitoring, 2-stage pressure reduction, 2nd stage with redundant design, pneum. priority switching in the event of power failure, with pressure gauges and emergency supply inlet (NIST), fully assembled on mounting plate incl. cover; incl. connection option for the reserve supply