


Kolibri, Rail AIR 0 / 8 - NIST



Manufacturer: Greggersen

Product 902606

Number:

Description

Kolibri, AIR, 8 l/min fixed

RAIL-MOUNTED UNIT NIST

Medical Compact Flowmeter

according to DIN EN ISO 15002; Flow fixed
to 8 l/min; Gas inlet NIST integrated in
rail clamp;

Outlet 9/16-18 UNF and Tubing Nipple

Gas type: Compressed Air AIR

Connection: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF / Tubing Nipple

Flow rate: 8 l/min

Flow settings: 0 + 8 l/min

Properties

Gas Type	AIR Compressed Air
----------	--------------------
