

Flowmeter, DIN, O₂, 0-15l/min



Manufacturer: Greggersen

Product 900753

Number:

Description

Flowmeter, O₂, 0...15 l/min

- PLUG-IN UNIT -

Medical flowmeter, pressure-compensated according to
DIN EN ISO 15002, with extension, fine
adjustment valve for metering the flow,
for coupling into a terminal unit

Gas type: Oxygen - O₂

Connection: Probe according to DIN 13260 Part 2

Overall length: 86.5 mm (probe front part to centre of flow outlet)

Outlet: 9/16-18 UNF

Flow rate: 0...15 l/min continuously adjustable

Inlet pressure: 450 kPa +/- 50 kPa

Properties

Gas Type	O2-Oxygen
----------	-----------

Inlet Standard	DIN
----------------	-----