


Kolibri, Rail O2 NEO - NIST



Manufacturer: Greggersen

Product 902603

Number:

Description

Kolibri, O₂, 0...1 l/min

RAIL-MOUNTED UNIT NIST

Medical compact flowmeter

according to DIN EN ISO 15002; flow rate is set via
handwheel in 10 flow settings;

gas inlet NIST integrated into rail clamp;

outlet 9/16-18 UNF and tubing nipple

Gas type: Oxygen O₂

Connection: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF / tubing nipple

Flow rate: 0...1 l/min

Flow settings: 0 / 0.1/ 0.15/ 0.2/ 0.25/ 0.3/ 0.4/ 0.5/ 0.6/ 0.8/ 1.0 l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
