

Forano, Front Part, DIN, N2O, Colour-Neutral



Manufacturer: Greggersen

Product 905113

Number:

Description

Terminal Unit FORANO - Front Part DIN

Colour-Neutral - Version III from 08/2019

Gas type: Nitrous oxide / N2O

Labelling: nitrous oxide / N2O

Gas-specific socket assembly for mounting
on FORANO Base Part.

Socket assembly according to DIN 13260 Part 2
with parking and operating position.

Release via release ring.

Colour coding:

Colour-Neutral (black)

Release ring labelling: printed

Language: English

Cartridge incl. sealing elements

Technical data:

Gas type: Nitrous oxide / N₂O

Labelling: nitrous oxide / N₂O

Colour coding: Colour-Neutral (black)

Front diameter: Ø 42mm

Scope of delivery:

- Gas-specific socket assembly
- Release ring incl. printed labelling
- Cartridge with sealing elements

Properties

Gas Type	N ₂ O-Nitrous Oxide
-----------------	--------------------------------
