


Ventus 22 mm Valve Insert +NIST CO2



Manufacturer: Greggersen

Product 900795

Number:

Description

Valve insert for Ventus 22 mm

Gas type : CO2

Design complies with DIN EN ISO 7396-1 and HTM 02-01.

Valve assembly consisting of:

- Shut-off valve DN 20, G3/4"
- CU pipe 22mm inlet / outlet
- Pressure gauge 0...16 bar
- Pressure sensor -1...6
- NIST emergency inlet point upstream of the shut-off valve
- NIST emergency inlet point downstream of the shut-off valve
- Physical disconnection point
- Gas type labelling

Properties

Fulfilled Standard	DIN EN ISO 7396-1 & HTM 02-01
---------------------------	-------------------------------

Gas Type	CO2-Carbon Dioxide
-----------------	--------------------

NIST Emergency Inlet Point	2 NIST Emergency Supply Points
-----------------------------------	--------------------------------