

AGSS Forano AIR AP TYPE1L



Manufacturer: Greggersen

Product 902076

Number:

Description

AGSS FORANO Type 1L
for surface mounting

according to DIN EN ISO 9170-2, with ejector, exhaust air recirculation,
with visual indicator for operational monitoring, mounted on a
mounting plate, cover made of sheet steel, galvanised, powder-
coated (W: 130 x H: 175 x D: 80 mm)

Flow rate: min. 50 l/min at 500 kPa (5 bar)

Gas type: Compressed air / AIR

Inlet: Cu pipe Ø 8 mm

Outlet: Exhaust air Cu pipe Ø 15 mm

Properties

Gas Type AIR Compressed Air