


Pressure Regulator Unit, AIR, Size 2, 8bar



Manufacturer: Greggersen

Product 325818

Number:

Description

Pressure regulator unit for reducing the vessel pressure to 8 bar, design in accordance with DIN EN ISO 7396-1, each pressure regulator stage duplicated with shut-off valves, safety valve 11 bar, contact sensor, completely mounted on a mounting plate.

Size: 2

Capacity: 300 Nm³/h

Gas type: Compressed air 8 bar / AIR-800

Pressure: 800 kPa (8 bar)

Contact sensor: 600 / 1000 kPa (6 bar / 10 bar)

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------