

Hose, VAC, 1.5m, NF - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 902529

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe

Hose and device connection NIST

Specification:

Gas type: Vacuum / VAC

Inlet: Probe according to NF S 90-116

Hose length: 1.5 m

Colour coding: Colour-Neutral

Outlet: gas-specific screw connection NIST

Properties

Gas Type

VAC-Vacuum