


Falke with Flow Paediatric O2 DIN Tubing Nipple +1DIN



Manufacturer: Greggersen

Product 910206

Number:

Description

Falke with Flow, Paediatric
with 1 quick-release coupling,
compact pressure regulator with contents gauge;
stepwise adjustable flow rate,
DIN quick-release coupling at a pressure of 450 kPa.

Gas type: Oxygen / O₂

Connection: Cylinder connection for gas cylinders to DIN 477, No. 9

Outlet: Tubing nipple / 1 quick-release coupling to DIN 13260 Part 2

Flow rate: 0...5 l/min tubing nipple - 120 l/min quick-release coupling

Flow settings: 0/0.1/ 0.3/ 0.5/ 0.6/ 0.8/ 1/ 2/ 3/ 4/ 5 l/min

Properties

Gas Type O₂-Oxygen