


Hose, VAC, 5m, DIN - DIN Coupling Colour-Neutral



Manufacturer: Greggersen

Product 902533

Number:

Description

Low-Pressure Hose
to DIN EN ISO 5359

Angled Probe,
Hose and coupling to DIN

Specification:

Gas type: Vacuum / VAC

Inlet: Probe to DIN 13260 Part 2

Hose length: 5.0 m

Colour coding: Colour-Neutral

Outlet: Coupling to DIN 13260 Part 2

Properties

Gas Type VAC-Vacuum