

Hose, VAC, 1.5m, DIN - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 900039

Number:

Description

Low-Pressure Hose

acc. to DIN EN ISO 5359

Angled Probe,

Hose and NIST device connection

Specification:

Gas type: Vacuum / VAC

Inlet: Probe acc. to DIN 13260 Part 2

Hose length: 1.5 m

Colour coding: Colour-Neutral

Outlet: gas-specific NIST screw connection

Properties

Gas Type

VAC-Vacuum