

Hose, CO2, 5m, DIN - NIST, ISO



Manufacturer: Greggersen

Product 900469

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow limiter,
hose and device connection NIST

Specification:

Gas type: Carbon dioxide / CO2

Inlet: Probe according to DIN 13260 Part 2

Hose length: 5.0 m

Colour coding: ISO 32 (grey)

Outlet: gas-specific NIST screw connection

Properties

Gas Type CO2-Carbon Dioxide