


AGSS FORANO AIR DVE2 TYPE1L



Manufacturer: Greggersen

Product 902077

Number:

Description

AGSS FORANO Type 1L

for ceiling pendants with integrated ejector

according to DIN EN ISO 9170-2, with exhaust air recirculation,
visual indicator for operational check and with
mounting bracket and outlet pipe.

Flow Rate: min. 50 l/min at 500 kPa (5 bar)

Gas Type: Compressed Air / AIR

Inlet: Tubing Nipple Ø 7 mm

Outlet: Exhaust Air Tubing Nipple Ø 15 mm

Properties

Gas Type

AIR Compressed Air