Compressed-Air Vacuum Regulator Specht -90 Probe ANTI



Manufacturer:	Greggersen
Product Number:	500244

Description

Vacuum Regulator Specht-90, DIN Plug-in Unit, non-magnetic

Medical compressed air suction unit for vacuum generation,
with Vacuum Gauge, fine dosing and quick-action shut-off valve,
silencer and filter holder.
Technical data:
Gas type: Compressed air / AIR
Inlet: Probe according to DIN 13260 Part 2
Control range: 0 to -85 kPa (0 to -0.85 bar)
Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH