Falke with Flow Neo O2 AGA Tubing Nipple



Manufacturer:	Greggersen
Product Number:	910336

Description

Falke with Flow, NEO,

Compact pressure regulator with high-pressure cylinder connection

acc. to SS 367615, contents gauge, incrementally adjustable

flow rate.

Gas type: Oxygen / O2

Connection: High-pressure cylinder connection acc. to SS 367615; No. E1 (W21.8)

Outlet: Tubing nipple at bottom

Flow rate: 0...1 l/min.

Flow settings: 0/ 0.1/ 0.15/ 0.2/ 0.25/ 0.3/ 0.4/ 0.5/ 0.6/ 0.8/ 1.0 l/min

Properties

Gas Type	O2-Oxygen

© 2026, Greggersen Gasetechnik GmbH