
Angled Probe AIR Hose Barb 7mm, DIN,ISO V4A



Manufacturer:	Greggersen
Product Number:	505151

Description

Angled Probe for Compressed Air

with 7 mm hose barb and non-return valve

Probe front part made of 1.4404 (stainless steel)

Gas type: Compressed Air / AIR

Connection: Probe according to DIN 13260 Part 2

Outlet: 7 mm Hose Barb

Colour coding: according to ISO 9170-1 (black/white)

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH