Vacuum Regulator Specht-90, AIR without Probe



Manufacturer:	Greggersen
Product Number:	500258

Description

Compressed-Air Vacuum Regulator Specht-90, connection G 1/4" male

Medical compressed-air vacuum regulator for vacuum generation,
with vacuum gauge, fine adjustment and quick-action shut-off valve,
silencer and filter holder.
Technical data:
Gas type: Compressed air / AIR
Inlet: G 1/4" male
Control range: 0 to -85 kPa (0 to -0.85 bar)
Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH