Vacuum Reg. Skua 2.0 -90, AIR,PROBE,DIN



Manufacturer:	Greggersen
Product Number:	904600

Description

Compressed-Air Vacuum Regulator Skua -90, Plug-in Unit DIN

Medical compressed-air vacuum regulator for dosing vacuum.
The vacuum is generated by means of compressed air according to the Venturi principle.
With 360° rotatable vacuum gauge; large, ergonomic
adjustment handwheel; rotatable exhaust sleeve.
The adjustment mechanism enables stepless, precise and
intuitive vacuum regulation. The rotatable
vacuum gauge allows the user to read the
set value from an optimal viewing angle. The
exhaust sleeve also allows the airflow to be individually aligned.
The regulator unit meets the requirements of
DIN EN ISO 10079-3 (high vacuum / high flow).
Material:
Housing: Aluminium, anodised
Handwheel: Plastic
Probe Front Part: Stainless steel
Technical Data:

© 2026, Greggersen Gasetechnik GmbH

Gas type: Med. compressed air / AIR

Inlet: Probe according to DIN 13260 Part 2

Control range: 0 to -85 kPa ± 10%

Flow rate: >20 l Freeflow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH