# Vacuum Reg. Pirol 2.0 -90,VAC,RAIL,NIST




| | |
|---|---|
| **Manufacturer:** | Greggersen |
| **Product Number:** | 904605 |

## Description

Vacuum Regulator Pirol -90, Rail-Mounted Unit NIST

Medical vacuum regulator for dosing vacuum.
With 360° rotatable vacuum gauge; large, ergonomic
adjustment handwheel. The adjustment mechanism enables stepless,
precise and intuitive regulation of the vacuum. The rotatable
vacuum gauge allows the user to read the
set value from an optimal viewing angle.
The regulator unit meets the requirements of
DIN EN ISO 10079-3 (high vacuum /high flow).
Material:
Housing: Aluminium, anodised
Handwheel: Plastic
Rail clamp: Aluminium, anodised
Technical data:

© 2026, Greggersen Gasetechnik GmbH

Gas type: Vacuum / VAC

Inlet: NIST housing according to DIN EN ISO 18082

Control range: 0 to -90 kPa ± 10%

Flow rate: > 20 l Freeflow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

## Properties

| | |
|---|---|
| **Gas Type** | VAC-Vacuum |

© 2026, Greggersen Gasetechnik GmbH