Secretion Overflow Protection 0.25 l VAC



Manufacturer:	Greggersen
Product Number:	900913

Description

Secretion Overflow Safety Device VAC 0.25 L
suitable for Skua + Pirol and Spatz

Connection: 9/16-18 UNF union nut
black, PA plastic
Outlet: Ø 6 - 10mm hose barb
Capacity: 0.25 L - APEC plastic
0.25 l complete with overflow safety device
and hose barb suitable for all Pirol, Skua and
Spatz regulators, connection from above.
Overflow safety device made of brass, chrome-plated

© 2026, Greggersen Gasetechnik GmbH