Vacuum Regulator Spatz -90 Probe DIN ANTIM



Manufacturer:	Greggersen
Product Number:	500243

Description

Vacuum Regulator Spatz-90, Plug-in Unit DIN, non-magnetic

Medical vacuum regulator based on the bypass principle,
with vacuum gauge, fine-adjustment and quick-shut-off valve.
Technical data:
Gas type: Vacuum / VAC
Inlet: Probe to DIN 13260 Part 2
Control range: -20 to -85 kPa (-0.2 to -0.85 bar)
Outlet: 9/16-18 UNF with vacuum tubing nipple

© 2026, Greggersen Gasetechnik GmbH