
Forano, Front Part, AGA, O2, ISO



Manufacturer:	Greggersen
Product Number:	905140

Description

Terminal Unit FORANO - Front Part AGA

Gas type: Oxygen / O2

Labelling: oxygen / O2

Gas-specific socket assembly for mounting on base part FORANO. Design complies with DIN EN ISO 9170-1

Socket assembly according to SS 87 524 30

Release via release ring.

Colour coding according to ISO 32.

Release ring labelling: printed

Cartridge incl. sealing elements

Technical data:

© 2026, Greggersen Gasetechnik GmbH

Gas type: Oxygen / O2

Labelling: oxygen / O2

Colour coding: ISO 32

Front diameter: Ø 42mm

Scope of delivery:

- gas-specific socket assembly
- release ring incl. printed text
- cartridge with sealing elements

Properties

Gas Type	O2-Oxygen

© 2026, Greggersen Gasetechnik GmbH