# Hose, AIR, 1.5m, DIN - NIST, Colour-Neutral



| | |
|---|---|
| **Manufacturer:** | Greggersen |
| **Product Number:** | 900040 |

## Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow limiter,

hose and NIST equipment connection

Specification:

Gas type: Compressed Air / AIR

Inlet: Probe according to DIN 13260 Part 2

Hose length: 1.5 m

Colour coding: Colour-Neutral

Outlet: gas-specific NIST screw connection

## Properties

| | |
|---|---|
| **Gas Type** | AIR Compressed Air |

© 2026, Greggersen Gasetechnik GmbH