Secretion Overflow Trap 0.25 l AIR



Manufacturer:	Greggersen
Product Number:	900911

Description

Secretion Overflow Safety Device AIR 0.25 L

Connection: 9/16-18 UNF Union Nut
Outlet: Tubing Nipple
Capacity: 0.25 l - APEC plastic
0.25 l compl. with overflow safety device
and tubing nipple suitable for all compressed-air-
operated regulators (SPECHT,
side connection to regulator)

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH