Hose, AIR, 2m DIN-NIST, Colour-Neutral angled



Manufacturer:	Greggersen
Product Number:	505097

Description

Special Low-Pressure Hose

++ MukoStar ++

Hose for supplying the device
MukoStar with compressed air 5 bar from a terminal unit;

Angled Probe with hose
with Non-Return Valve, hose diameter 6.7mm
and device connection NIST - 90° angled

Specification:
Gas type: Medical Air / AIR
Inlet: Probe according to DIN 13260 Part 2
Hose length: 2 m
Colour coding: Colour-Neutral
Outlet: gas-specific screw connection NIST

© 2026, Greggersen Gasetechnik GmbH

90° angled

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH