Vacuum Reg. Pirol 2.0 -90,VAC,PROBE,DIN



Manufacturer:	Greggersen
Product Number:	904604

Description

Vacuum Regulator Pirol -90, Plug-in Unit DIN

Medical vacuum regulator for vacuum dosing.
With 360° rotatable vacuum gauge; large, ergonomic
adjustment handwheel. The adjustment mechanism enables stepless,
precise and intuitive vacuum regulation. The rotatable
vacuum gauge allows the user to read the
set value from an optimal viewing angle.
The regulator unit meets the requirements of
DIN EN ISO 10079-3 (high vacuum / high flow).
Material:
Housing: aluminium, anodised
Handwheel: plastic
Probe Front Part: stainless steel
Technical Data:
Gas type: Vacuum / VAC
Inlet: Probe according to DIN 13260 Part 2
Control range: 0 to -90 kPa ± 10%

© 2026, Greggersen Gasetechnik GmbH

Flow rate: > 20 l free flow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	VAC-Vacuum

© 2026, Greggersen Gasetechnik GmbH