Nebuliser Plug-in Valve with Clamp AIR



Manufacturer:	Greggersen
Product Number:	500370

Description

Fitis Nebuliser Plug-in Valve AIR

Plug-in Unit

Preset nebuliser plug-in valve, with side clamp,

for coupling into a terminal unit,

with Tubing Nipple

Technical Data:

Gas type: Compressed Air / AIR

Connection: Probe to DIN 13260 Part 2

Outlet: 9/16-18 UNF/Tubing Nipple

Flow Rate: 5 l/min ± 20%

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH