Vacuum Reg. Pirol 2.0 -30,VAC,RAIL,CH



Manufacturer:	Greggersen
Product Number:	505656

Description

Vacuum Regulator Pirol -30, Rail-Mounted Unit CH, NIST
with spring-loaded rail clamp for up to 50 mm rails

Medical vacuum regulator for dosing vacuum.
With 360° rotatable vacuum gauge; large, ergonomic adjustment handwheel; rail-mounted unit incl. spring-loaded rail clamp.
The gas inlet is integrated into the rail clamp (NIST).
The adjustment mechanism enables stepless, precise and intuitive vacuum regulation. The rotatable vacuum gauge enables the user to read the set value from an optimum viewing angle.
The regulator unit meets the requirements of DIN EN ISO 10079-3 (high vacuum / high flow).
Material:
Housing: Aluminium, anodised

© 2026, Greggersen Gasetechnik GmbH

Handwheel: Plastic

Rail clamp: Aluminium, anodised

Technical data:

Gas type: Vacuum / VAC

Inlet: NIST housing according to DIN EN ISO 18082

Control range: 0 to -30 kPa

Flow rate: >20 l free flow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

© 2026, Greggersen Gasetechnik GmbH